



Data Collection and Preprocessing Phase

Date	15 March 2024
Team ID	SWTID1720425899
Project Title	Covid Vision: Advanced COVID-19 Detection from Lung X-rays with Deep Learning
Maximum Marks	6 Marks

Preprocessing:

Section	Description	
Data Overview	The dataset contains data of chest X-ray images for Covid- 19 positive cases along with Normal and Viral Pneumonia images.	
Resizing	Resize images to a target size of (256, 256)	
Normalization	rescale=1./255: This parameter is used to normalize pixel values of images. In many image datasets, pixel values range from 0 to 255.	
Data Augmentation	Applied augmentation techniques such as zooming range = 0.2, shearing range = 0.2.	
Data Preprocessing Code Screenshots		
Loading Data	<pre>!pip install opendatasets import opendatasets as od import pandas as pd od.download("https://www.kaggle.com/datasets/tawsifurrahman/covid19-radiography-database")</pre>	





	▶ [plp install opendatasets inport opendatasets as od inport pandatasets as od inport pandas as pd od. dom.lond/[https://www.kagle.com/datasets/taxelfurralman/covid19-radiography-database*) ➡ Requirement already staisfied: constitution as a part of the part o
Resizing	IMAGE_SIZE = (256, 256) IMAGE_SHAPE = IMAGE_SIZE + (3,)
Normalization	<pre>from tensorflow.keras.preprocessing.image import ImageDataGenerator img_height, img_wldth= IMAGE_SIZE batch_size=16 train_datagen = ImageDataGenerator(rescale=1./255,zoom_range=0.2,shear_range=0.2,validation_split=0.3)</pre>
Data Augmentation	[] from tensorflow.keras.preprocessing.image import ImageDataGenerator img_height, img_width= IMAGE_SIZE batch_size=16 train_datagen = ImageDataGenerator(rescale=1./255,zoom_range=0.2,shear_range=0.2,validation_split=0.3)